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Bone and collagen matrix OsteoBiol® English version -

Class III medical device, bone and collagen matrix of heterologous origin, for the replacement of the autologous bone represented by the matrix and the porous form.

Putty, bone paste and collagen of heterologous origin. The collagen component consists of a matrix of type I and type III collagen and polyunsaturated fatty acids. PARTICLE SIZE: <0,300 mm; the product is slightly radioopaque. PACKAGING: sterile disposable syringes.

CLINICAL INDICATIONS: filler for non-infectious, non-sclerotic, vascularized bone defects, on non-load bearing sites.

DIRECTIONS FOR USE: apply the necessary amount of product directly to the receiving site. OsteoBiol® Putty must be grafted only in cavities that can ensure its stability. Therefore, it must not be grafted in two wall defects or in lateral access sinus lift procedures. It is possible to mix the product with bone substitutes in granules, as needed. In Dentistry, it is advisable to protect the graft with a resorbable barrier membrane to increase stability before proceeding with flap surgery. Note for filling large cavities with volume greater than 1 cm³, mix one part of product with three parts of bone substitutes in granules (final volume of Putty: 25%).

SIDE EFFECTS: none known. Although the grafting procedure adopted has demonstrated the inhibition of adverse antibody reactions, the remote possibility of graft rejection cannot be excluded.

CONTRADICTIONS: do not use in patients who have pathologies or conditions that contraindicate surgery or who have an acute or chronic infection (osteomyelitis). Not indicated for patients with immunodeficiency or suffering from autoimmune diseases, particularly children and the elderly. Do not use directly or near large vessels or nerves. Do not use on poorly vascularized, necrotic or infected receiving sites. Do not use on sites that do not guarantee the stability and protection of the product after grafting. The product must not be used if the packaging is damaged or opened for use.

WARNING: This product has not been tested on pregnant women. The product is for single use only and, as such, should not be sterilized. The graft site must be properly washed and cleaned, and the product must be well stabilized and protected.

Product reserved for surgeons with experience in biomaterial implants; strictly adhere to the instructions given in the scientific literature on heterologous grafts, especially the scientific literature on OsteoBiol® devices available on the website www.osteobiol.com, as well as the instructions hereto. For every single application, the specialist surgeon has the responsibility to ensure proper clinical management of the biomaterials and the patient.

EXPIRY DATE: the product is valid for five years from the date of sterilization and the expiry date is indicated on the label.

STORAGE: The product must be stored away from direct heat sources and sunlight, in a dry and clean place, under the conditions indicated on the label.

DISPOSAL: the product does not contain hazardous substances and it may be disassembled and disposed of as regular municipal waste.

PRODUCT ORIGIN: the last letter on the label of the product code REF indicates the source of the material (S-porcino; E-equina).

Symbols used on the label:

① single use ② year and month of expiry ③ lot number ④ carefully read the instructions ⑤ Sterile. Sterilization performed by irradiation with gamma rays. If the packaging is open or damaged the sterility of the product is no longer guaranteed.

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DIRECTIONS FOR USE: apply the necessary amount of product directly to the receiving site. OsteoBiol® Putty must be grafted only in cavities that can ensure its stability. Therefore, it must not be grafted in two wall defects or in lateral access sinus lift procedures. It is possible to mix the product with bone substitutes in granules, as needed. In Dentistry, it is advisable to protect the graft with a resorbable barrier membrane to increase stability before proceeding with flap surgery. Note for filling large cavities with volume greater than 1 cm³, mix one part of product with three parts of bone substitutes in granules (final volume of Putty: 25%).

SIDE EFFECTS: none known. Although the grafting procedure adopted has demonstrated the inhibition of adverse antibody reactions, the remote possibility of graft rejection cannot be excluded.

CONTRADICTIONS: do not use in patients who have pathologies or conditions that contraindicate surgery or who have an acute or chronic infection (osteomielitis). Not indicated for patients with immunodeficiency or suffering from autoimmune diseases, particularly children and the elderly. Do not use directly or near large vessels or nerves. Do not use on poorly vascularized, necrotic or infected receiving sites. Do not use on sites that do not guarantee the stability and protection of the product after grafting. The product must not be used if the packaging is damaged or opened for use.

WARNING: This product has not been tested on pregnant women. The product is for single use only and, as such, should not be sterilized. The graft site must be properly washed and cleaned, and the product must be well stabilized and protected.

Product reserved for surgeons with experience in biomaterial implants; strictly adhere to the instructions given in the scientific literature on heterologous grafts, especially the scientific literature on OsteoBiol® devices available on the website www.osteobiol.com, as well as the instructions hereto.

EXPIRY DATE: the product is valid for five years from the date of sterilization and the expiry date is indicated on the label.

STORAGE: The product must be stored away from direct heat sources and sunlight, in a dry and clean place, under the conditions indicated on the label.

DISPOSAL: the product does not contain hazardous substances and it may be disassembled and disposed of as regular municipal waste.

PRODUCT ORIGIN: the last letter on the label of the product code REF indicates the source of the material (S-porcino; E-equina).

Symbols used on the label:

① single use ② year and month of expiry ③ lot number ④ carefully read the instructions ⑤ Sterile. Sterilization performed by irradiation with gamma rays. If the packaging is open or damaged the sterility of the product is no longer guaranteed.

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