**AB GERMA®**

## **Vacuum Positioning Pillow**

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# Directions for use

## 2. Precautions

- On a current basis check that the vacuum pump works correctly and that it is possible to get the air out of the Pillow to make it hard.
- Avoid handling the pillow next to sharp objects. If the pillow is pulled be aware of obstacles, sharp edges etc.
- Avoid handling the pillow next to open fire.
- On a current basis check that the pillow has no cracks.
- Own modifications on products not allowed – Your medical specialist to be responsible for any modification.
- When decided to take the product out of service and destroy – make sure to follow your local instructions or contact us for additional information.
- DentaComfort pillow made in compliance with MDR 2017/745 Class 1 concerning Medical Device.

### Introduction Vacuum Pillow

- How to use your Germa Vacuum pillow.

The Gernia Vacuum pillow is intended to be used as a tool for, positioning and fixation. Immobilization or stabilisation of patients in many different situations. In order to achieve the desired form and rigidity of the pillow it is necessary to mold it and gradually create vacuum in the pillow by evacuating the air. This can be done by using Germa Vacuum Pump or another kind of vacuum source. The pillow consists of an outer shell of airtight material and a textile inner bag containing granules.

### Precautions and test

On a current basis check that the vacuum pump works correctly and that it is possible to evacuate the air from the pillow. Avoid handling of the vacuum pillow on or close to sharp objects or next to open fire. Check for damages, cracks or wear and tear.

### Specifications

The vacuum pillow tested and approved according to EN 1021-1. Materials are cover material in PVC-film or PVC coated fabric, inner bag in polyamide and styropor granules. Sizes are: 80 x 100 cm, 80 x 200 cm, 55 x 30 cm, 55 x 65 cm, 55 x 85 cm and 40 x 50 cm and optionally made to customers measurements. Recommended temperature of use -20° C - +40° .

### Storage

Recommended storage temperature -20 ° C - +70 °.

### In use

Open the valve and put the unit on an even surface. When you have obtained an even distribution of the granules. Shape and fold the unit in its initial position.

Put the vacuum pillow in the desired position on the patient. For pillows with t-valve: Open the valve, evacuate the air and close the valve. For screw valve: Close the valve until it is totally closed – anti clock. If you are not satisfied with the support and stabilization of the patient. Just open the valve, relocate the pillow in desired position, close the valve, and evacuate again.

### Cleaning and disinfecting

After use of the product it can be necessary to rinse it in soap and water. Allow the pillow to dry completely before storing it. If necessary to disinfect. wipe the complete pillow with a solution containing of 70 % alcohol or other most common disinfectants e.g. Lysetol AF.

### Maintenance and repair

Test the pillow and pump at regular intervals to confirm the necessary airtightness and functions. Eventual leaks can be found by visual inspection or by reduced filling of air into the pillow followed by

covering it with soapy water all smear the surface and the welds watching out for air bubbles from the pillow. If leaks have been found, clean the area and let it dry. Find a patch in a suitable size and put glue on the patch as well as on the leakage. Wait approx. 5 minutes and place the patch on the leakage. Firmly press the whole area of the patch establishing good contact and press all air out under the patch. The repair will receive full strength after 24 hours.

### **Test**

For testing evacuate the air and create rigidity, close the valve and leave the pillow for about two hours. If it is still as rigid as it was when you evacuated, the pillow may be considered tight. If the pump is not evacuating the pillow, check the pump by closing the tube connector with your fingertip and pump a stroke two-three times. A negative pressure should be created in the tubing and can be felt at your fingertip. Repair kits for pillow is available.

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### **Discard of product**

DentaComfort cushions can be disposed of with normal household waste as specified by local recycling recommendations for plastics according to their active duration. Materials used in the packaging are environmentally friendly and can be recycled.

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**Report no. :**

(To be assigned in by  
Manufacturer)

**Please print or type or fill in by MS Word, all sections a. – i. Fax or e-mail to your supplier.  
Please use repair order when reporting repairs not associated with complaints.**

For claims please specify :  
Purchase Date

**Description of product**

**Description of Incident/Problem:** 1) What, 2) How?

When? In connection with: 1) transport, 2) unpacking, 3) preparing, 4) test, 5) cleaning, 6) normal use, 7) training, 8) disassembly, 9) assembly, 10) storage, 11) other:

**Info on the product, prior to problem,**

Storage: 1) stock, 2) site of use, 3) packed, 4) unpacked, 5) softpack, 6) other:

Used cleaning and disinfecting agents/Sterilizing method:

Pattern of use: 1) stand-by for emergency, 2) routine use, 3) stand by, 4) other:

Environment: 1) Indoor, 2) Outdoor, 3) Heat, 4) Cold, 5) Humidity, 6) Sunlight, 7) Gases, 8) other:

---

Which **actions** have been: a) taken, b) planned by the reporter or user at the site of incident?

1) examination of the product on site, 2) replacement of the product, 3) other:

---

Suggested actions by the distributor/complainant

---

---

**Returning of samples/faulty product(s),**

- 1) Already shipped/Will be shipped to \_\_\_\_\_,
- 2) Pending request from manufacturer,
- 3) Product/parts not available:

---

Date:

---

**We acknowledge the receipt of this report by manufacturer and will start processing it immediately. Please use the assigned complaint report number for further reference.**

**Your contact person in this matter is:** \_\_\_\_\_

Name:

Direct Tel.

Direct Fax.

Init/Date:

---

# DentaComfort pillow

The DentaComfort pillow is used by dentists for the positioning of patients in the dental chair.

**DentaComfort** vacuum pillow adjusts to the patient head shape on the top and the chair on the back and helps patient to relax and lie still in desired position.

The DentaComfort pillow has a range of benefits and features for both dentists and patients.

**Dentist ref:**

- My patient was very calm and relaxed throughout the treatment .

Easy to give the patient comfortable support in a variety of positions, quickly

Face of the pillow shaped after patient head.


The back of the pillow locking position.

Hyperextension effortless for the patient

**Parent:**

- Our children feel more secure and enjoy going to the dentist now.

Use your saliva suction for any second stabilizing of the cushion and it is just as quick to reposition.

**Patient ref:**

- Just fantastic support to my neck, very relaxing and safe.

**GERMA**

AB GERMA  
Box 2127, SE-291 02 KRISTIANSTAD, SWEDEN  
Phone: +46 (0)44 12 30 30

Fax: +46-(0)44 10 31 79  
E-mail: info@germa.se  
Internet: www.germa.se

## DentaComfort Pillow

**Manufacturer:**

AB Germa  
Industrigatan 54-56  
SE-291 36 Kristianstad  
Sweden

Tel: +46 (0)44 12 30 30  
Fax: +46 (0)44 10 31 79

E-mail: [info@germa.se](mailto:info@germa.se)

[www.germaprotec.com](http://www.germaprotec.com)

### Directions for use

#### Vacuum DentaComfort Pillow

The DentaComfort pillow equipment for dental practices is designed to help, support and stabilize the patients head and neck during treatment. It is designed to be used on most medical dentist chairs and offer ergonomic workplace during complex treatments and offer additional safety plus maximum of relaxing and comfort for the patient.

To create desire shape and support with DentaComfort pillow place the pillow under the patient´s head and neck and create the desired position, create vacuum inside the pillow by using a saliva sucker to make the pillow form stable in the desired position.

#### Storage

Store the pillow in room temperature and dry place.

#### Specification

DentaComfort pillow is made of PVC film, latex free outer cover with slide valve. Detached inner bag of Nylon filled with granulates of polystyrene. Inner bag is designed to control granulates and to prevent any loss of granulates if there arise a damage on the outer material.

#### Safe usage

1 Place the pillow under the patient´s head and neck and create the desired shape and position by lifting gently on the side with both hands.


2 Open the slide valve and suck out air with the saliva sucker until the desired stability is achieved by the vacuum. It is very quick to take out air.

Do not push too hard on valve towards the pillow, this to make sure you have free air flow.

3 Make sure to close the slide valve by sliding the valve to lock position and remove the saliva sucker. You have created the desired position; note hardness of the pillow is determent by how much air you take out from the pillow.

### **Warnings and precautions**

**Important**, do not place the pillow in front of or on the face of a patient.

Make sure not to have the pillow close to any sharp tools, edges objects, strong heat or fire. Always avoid damaging the outside cover.

Inspect the pillow for any possible damages before use and make sure the pillow and valve is in good condition for use. Do not use the pillow if you have discovered any damages.

### **Cleaning and disinfection**

Clean the pillow with water and soap after use. When called for disinfection it can be cleaned with most common disinfectants. Dry the pillow before next usage.

Operational room covers shall be washed according to tags and disinfections regulations.

### **Warranty**

Warranty 1 (one) year from delivery date. No liability for outer violence or damages from improper handling.

### **Discard of product**

DentaComfort pillows may be scraped with regular household garbage at the end of its duration. Material used in packing is environmentally friendly and can be recycled.

### **Accessories**

Covers, operational room suitable cotton covers and single use cover (will be introduced)