UNIVET

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(EN) Instructions for use - MAGNIFYING SYSTEM:

- TTL BINOCULAR LOUPES
- FLIP UP BINOCULAR LOUPES
- FLIP UP BINOCULAR HEADGEAR

Thank you for your trust in us and congratulations on your choice.

Read the following information before using the device.

DESCRIPTION OF THE COMPONENTS

The user guide for the magnifying system, the list of components and spare parts are indicated in the documents provided along with the device.

SYMBOLS

In compliance with Regulation 2017/745 on medical devices

Date of manufacture Manufacturer Catalogue number

Serial number Batch number Instructions for use Caution

INTENDED USE

The device is intended to improve visual acuity in a narrow range of the visual field in medical applications and intended to compensate for a visual defect when equipped with ophthalmic load-bearing lenses.

Intended users: medical professionals.

Intended patients: no limitation except as stated in the "contraindications" section.

DEVICES THAT CAN BE USED TOGETHER

The device can be comfortably attached to all Univet Lighting system models, thus providing an integrated vision system.

SPECIFICATIONS

TYPE OF MAGNIFYING SYSTEM	TYPE OF OPTICAL SYSTEM	MAGNIFICATION	WORKING DISTANCE [mm]
TTL BINOCULAR LOUPES	GALILEAN	2.0X PRO HD	300
		2.5X PRO HD	350
		3.0X PRO HD	400
		3.5X PRO HD	450
		3.5X PRO HD	500
TTL BINOCULAR LOUPES	PRISMATIC	3.5X XS ULTRA-HD	300
		4.0X XS ULTRA-HD	350
		4.0X XS ULTRA-HD	400
		4.5X XS ULTRA-HD	450
		5.0X XS ULTRA-HD	500
FLIP UP BINOCULAR LOUPES	GALILEAN	2.5X PRO HD	300
			350
			400
	PRISMATIC	3.5X ULTRA-HD	300
		4.5X ULTRA-HD	350
FLIP UP BINOCULAR HEADGEAR	PRISMATIC	6.0X ULTRA-HD	300
			350
			350
			400
			400
FLIP UP BINOCULAR HEADGEAR	PRISMATIC	3.5X ULTRA-HD	300
		4.5X ULTRA-HD	350
		4.5X ULTRA-HD	400
		4.5X ULTRA-HD	450
		6.0X ULTRA-HD	300
	350		
	400		

CONTRAINDICATIONS

The device (loupes equipped with TTL magnifying system) is not recommended for users suffering from strabismus.

FLIP UP magnifying systems are recommended for operators suffering from strabismus when

used along with the corrective spectacles thereof in use.

The device is contraindicated in the event of applications requiring image stability that cannot be achieved with devices mounted on the user's head.

WARNINGS

Read the following warnings carefully before using the device.

Make sure that it is suitable for the intended use.

Clean the device before use.

Not use the device if it has been damaged but contact your authorised dealer.

Use the supplied components only.

Ensure that the shields (if provided) are correctly fitted, otherwise stop the procedure in progress.

Do not tamper with or disassemble any component of the system.

Do not attempt any repairs.

In the event of a malfunction, stop using the device, contact your authorised dealer.

When not in use, store the device in its container.

PRECAUTIONS

For any component adjustment and fitting instructions when using for the first time, please see the information provided by the manufacturer.

Avoid contact with heat and cold sources or open flames.

To avoid damaging the eyes, please do not look directly into light sources.

Avoid sudden movements with your head when using the device.

Avoid walking with the device worn.

Do not wear the device over headgear.

Please use the elastic band when using loupes equipped with TTL and FLIP UP magnifying system.

ADVERSE EFFECTS

No adverse effects detected

SUPPLY MODE

This device is not supplied sterile. This is not a single use device.

STORAGE AND/OR HANDLING

Store the device at the temperatures and humidity indicated in the "Technical Specifications" section.

When not in use, store the cleaned device in its original packaging.

Handle the device strictly by holding by the temples (in the case of Loupes) or by the bands (in the case of Headgears).

MAINTENANCE

Routine maintenance consists of cleaning the device as described in the "cleaning instructions" section.

Ensure that the components of the device are fitted and adjusted correctly before use.

Keep the original packaging required for return/maintenance shipment and contact your authorised dealer, see the instructions laid down by the warranty.

TRAINING AND QUALIFICATIONS FOR USE

None.

CLEANING INSTRUCTIONS

Use a microfiber cloth to clean the lenses.

All the other surfaces can be cleaned with a soft cloth soaked with a gentle detergent solution.

This device cannot be sterilised. Do not clean in an ultrasonic bath.

DISPOSAL INSTRUCTIONS

The device does not contain toxic or hazardous materials. Disinfect the device before disposal.

In the event of disposal/demolishing, do not dispose of the device in the environment, comply with the relevant national regulations.

HEALTH RISKS

No health risk detected.

WARRANTY

See the instructions laid down by the warranty.

Report any serious accidents involving the device to the manufacturer and the competent authority of the Member state in which the user and/or the patient resides.

TECHNICAL SPECIFICATIONS

Operating conditions

Room temperature.....from +10°C to +40°C

Relative humidity.....from 10% to 90%

Storage and transportation conditions

Room temperature.....from -10°C to + 50°C

Relative humidity.....from 10% to 90%

Store in a dry place away from heat and cold

CLINICAL BENEFITS

Given that it offers greater visual acuity (magnification) and facilitates a smaller angle of the position of the operator's neck, the device helps the user to maintain the correct posture during use, hence avoiding the onset of musculoskeletal problems.

No adverse effects for the user.

Univet's magnifying systems are effective and safe when used as intended.

For languages not included in these use instructions, please see the following link:

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