SERALENE®**de GEBRAUCHSANWEISUNG****NAHTMATERIAL AUS POLYVINYLIDENFLUORID****en INSTRUCTIONS FOR USE****SUTURE MATERIALS OF POLYVINYLIDENE FLUORIDE****fr MODE D'EMPLOI****MATERIAU DE SUTURE EN FLUORURE DE POLYVINYLE****it ISTRUZIONI PER L'USO****MATERIALE DA SUTURA IN FLUORURO POLIVINILIDENICO****es INSTRUCCIONES DE USO****MATERIAL DE SUTURA DE FLUORURO DE POLIVINILIDENO****pt INSTRUÇÕES DE SERVIÇO****MATERIAL DE SUTURA À BASE DE FLUORETO DE POLIVINILO****nl GEBRUIKSAANWIJZING****HECHTMATERIAAL UIT POLYVINYLIDEENFLUORIDE****pl INSTRUKCJA OBSŁUGI****NICI CHIRURGICZNE Z POLIFLUORKU WINYLIDENU****cs NÁVOD K POUŽITÍ****CHIRURGICKÝ ŠICÍ MATERIÁL Z POLYVINYLIDENFLUORIDU****tr KULLANIM TALİMATI****POLİVİNİLİDEN FLORÜRDEN SÜTÜR MATERYALİ****lt VARTOJIMO INSTRUKCIJA****MEDŽIAGOS ŽAIZDOMS SIŪTI IŠ POLIVINILLIDENO FLUORIDO**

lv LIETOŠANAS INSTRUKCIJA

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SUTURE MATERIALS OF POLYVINYLIDENE FLUORIDE

- SERALENE®

WARNINGS AND PRECAUTIONS

- Read the instructions for use in their entirety before using the product.
- The product should only be used by qualified personnel, that is responsible to be familiar with both appropriate surgical techniques and anatomy prior to use.
- The product is delivered sterile (sterilization method: ethylene oxide).

- The product is intended for single use only.

- Do not resterilize or re-use the product as critical changes in performance and safety are possible.

- Do not use after expiration date.

- Do not use if the product shows signs of damage (e.g. cuts or kinks) or if the outer packaging or sterile packaging has been damaged or unintentionally opened prior to use.

- Use the product only for the intended indication.

DESCRIPTION AND INTENDED PURPOSE

SERALENE® is a synthetic monofilament non-absorbable surgical suture material. The basic material of the thread is polyvinylidene fluoride. SERALENE® is dyed with solvent blue, Colour Index No. 61568. SERALENE® is available in various gauge sizes and lengths.

It is supplied with various types of needles made of surgical steel or without needles. It may be supplied in combination with other components, such as a suture guiding tube. Pledgets may be enclosed with SERALENE® to reinforce the body's own structures.

The suture material complies with the current edition of the European Pharmacopoeia monograph (EP) "Sutures, sterile non-absorbable" (fila non resorbilia sterilia)".

SERALENE® is suitable for adaptation of soft tissues and for ligatures where non-absorbable suture material is indicated, including use in cardiovascular and neurological procedures as well as in microsurgery and ophthalmic procedures. Polyvinylidene fluoride threads can also be used as holding sutures and for marking purposes.

INSTRUCTIONS FOR USE

In the choice and use of suture material, account should be taken of the condition of the patient, the experience of the qualified personnel, the surgical technique, and the type

and size of the wound. Standard knot-tying techniques should be used to ensure knot-holding security. In order to prevent needle damage during suturing it is recommended that the needle be grasped at a point between the end of its first third (as seen from its end) and its midpoint.

Unneeded suture material may be used for ligation or may be connected with surgical eye needles. For suitable needles and application instructions, please refer to the Instructions for Use of Serag-Wiessner sterile eye needles. The suitability of other manufacturers' products has not been tested. When using other brands, their specific instructions must be observed.

ACTION

The function of a medical suture in wound care is to attach tissue to tissue or tissue to foreign material. Ligatures are a special type of suture used to close hollow organs.

After a slight initial tissue reaction, suture material of polyvinylidene fluoride is slowly and progressively encapsulated by connective tissue. Because of its monofilament thread structure, SERALENE® can also be used in tissue at risk of infection.

CONTRAINDICATIONS

The device must not be used in patients with allergies to any of the product components.

RESTRICTIONS AND LIMITATIONS OF USE / NOTES / INTERACTIONS

The requirements of the suture material and the techniques required vary with the indication. In the choice of suture material, account should be taken of the characteristics of the suture material concerned. Like all foreign bodies, the device can give rise to calculus formation when in pro-longed contact with salt-containing solutions (e.g. in the biliary or urinary tract). In the case of sutures in infected tissue the use of suitable ancillary measures is recommended. Surgical instruments should be handled with care in order to avoid damage to threads. Bending of needles or grasping them outside of the recommended area can result in needle breakage. Because of the risk of infection from needle pricks, needles should be handled with particular care.

ADVERSE REACTIONS

The following adverse reactions can occur in association with correct use of SERALENE®: localized irritation, inflammatory foreign-body reactions, suture fistulas or granulomas.

Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

STORAGE CONDITIONS

The storage conditions apply as indicated on the label.

DISPOSAL INFORMATION

Disposal of unused products or waste material with particular care due to the potential risk of contamination and/or injury. In all other cases observe the specific national regulations.

MEANING OF THE SYMBOLS

Symbols are explained at the end of the instructions for use.

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SERAG-WIESSNER GmbH & Co. KG