OsteoBiol®

by Tecnos

Lamina

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Bone and collagen matrix OsteoBiol® - English version -

Class III medical device, bone and collagen matrix of heterologous origin, for the replacement of the autologous bone represented by the matrix and the porous form.

Lamina: preformed bone lamina-shaped heterologous bone matrix and collagen, osteoconductive, resorbable, biocompatible, hygroscopic. CLINICAL INDICATIONS: coadjutant for the reconstruction or the recuperation of total or partial osseous perdices, as a filler or protective means for non-infected, non-sclerotic and well vascularized bone defects. PACKAGING: sterile disposable blister.

DIRECTIONS FOR USE:

•LS10HSL51010E codes "Lamina Curva" (bent lamina) - do not hydrate before use: shape and adapt the product to the anatomical characteristics of the receiving site with scissors or sterile instruments and then fix it with titanium micro screws directly by suturing to the surrounding tissues with an atraumatic needle with triangular cross-section.

•Codes LS3L5L53L53LE : medium degree of mineralization with rigidity of the preformed good to the natural bone; shape the product by hydration in sterile conditions and fix it to the receiving site with synthetic instruments; then suture the flaps.

•All the other codes: decreased degree of mineralization, flexible preformed; open the package and hydrate the product by soaking it for 5'10 minute in sterile physiological solution. Once the product has acquired the desired plasticity, it can be placed on the site previously prepared and fixed with synthetic instruments; then you can suture the flaps.

SIDE EFFECTS: none known. Although the grafting procedure adopted has demonstrated the inhibition of adverse antibody reactions, the remote possibility of graft rejection cannot be excluded. CONTRAINDICATIONS: do not use in patients who have pathologies or conditons that contraindicate surgery or who have an acute or chronic infection (osteomyelitis). Not indicated for patients with immunodeficiency or suffering from autoimmune diseases, particularly children and the elderly. Do not use directly or near large vessels or nerves. Do not use on poorly vascularized, necrotic or infected receiving sites. Do not use on sites that do not guarantee the stability and protection of the product after grafting. The product must not be used if the packaging is damaged or opened before use. WARNING: This product has not been tested on pregnant women. The product is for single use only and, as such, should not be resterilized. The graft site must be properly washed and cleaned, and the product must be well stabilized and protected. Product reserved for surgeons with experience in biomaterial implants: strictly adhere to the instructions given in the scientific literature on heterologous grafts, especially the scientific literature on OsteoBiol® devices available on the website www.osteobiol.com, as well as the instructions herein. It is the responsibility of the specialist surgeon who has the responsibility to ensure proper clinical management of the biomaterials used in the patient. EXPIRY DATE: the product is valid for five years from the date of sterilization and the expiry date is indicated on the label. STORAGE: The product must be stored away from direct heat sources and sunlight, in a dry and clean place, under the conditions indicated on the label. DISPOSAL: the product does not contain hazardous substances and it may be disseminated and disposed of as regular municipal waste. PRODUCT ORIGIN: the last letter on the label of the product code REF indicates the source of the material (S = porcine, E = equine).

Simbols used on the label:

single use year and month of expiry lot number carefully read the instructions Steril. Sterilization performed by irradiation with gamma rays. If the packaging is open or damaged the sterility of the product is no longer guaranteed limitation of the temperature

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Bone and collagen matrix TecnoSS® - English version -

Class III medical device, bone and collagen matrix of heterologous origin, for the replacement of the autologous bone represented by the matrix and the porous form. **Lamina:** preformed bone lamina-shaped heterologous bone matrix and collagen, osteoconductive, resorbable, biocompatible, hygroscopic. **CLINICAL INDICATIONS:** an adjuvant for the reconstruction or partial or total recovery of lost bone portions, as a filler or protective means for non-infected, non-sclerotic and well vascularized bone defects. **PACKAGING:** sterile disposable blister. **DIRECTIONS FOR USE:** -LS10HS/LS10HE codes "Lamina Curva" (bent lamina) - do not hydrate before use: shape and adapt the product to the anatomical characteristics of the receiving site with scissors or sterile instruments and then fix it with titanium micro screws or directly by suturing to the surrounding tissues with an atraumatic needle with triangular cross-section. -Codes LS35LS/LS35LE/LS15LS/LS15LE/LN32LS/LN35HE: medium degree of mineralization with rigidity of the preformed close to the natural bone: shape the product by hydration in sterile conditions and fix it to the receiving site with synthetic instruments; then suture the flaps. -All the other codes: decreased degree of mineralization, flexible preformed; open the package and hydrate the product by soaking it for 5/10 minute in sterile physiological solution. Once the product has acquired the desired plasticity, it can be placed on the site previously prepared and fixed by means of synthetic instruments; then you can suture the flaps. **SIDE EFFECTS:** none known. Although the grafting procedure adopted has demonstrated the inhibition of adverse antibody reactions, the remote possibility of graft rejection cannot be excluded. **CONTRAINDICATIONS:** do not use in patients who have pathologies or conditions that contraindicate surgery or who have an acute or chronic infection (osteomyelitis). Not indicated for patients with immunodeficiency or suffering from autoimmune diseases, particularly children and the elderly. Do not use directly or near large

vessels or nerves. Do not use on poorly vascularized, necrotic or infected receiving sites. Do not use on sites that do not guarantee the stability and protection of the product after grafting. The product must not be used if the packaging is damaged or opened before use. **WARNING:** This product has not been tested on pregnant women. The product is for single use only and, as such, should not be re-sterilized. The graft site must be properly washed and cleaned, and the product must be well stabilized and protected. Product reserved for surgeons with experience in biomaterial implants: strictly adhere to the instructions given in the scientific literature on heterologous grafts, especially the scientific literature on TecnoSS® devices available on the website www.TecnoSS.com, as well as the instructions herein. For every single application, the specialist surgeon has the responsibility to ensure proper clinical management of the biomaterials and the patient. **EXPIRY DATE:** the product is valid for five years from the date of sterilization and the expiry date is indicated on the label. **STORAGE:** the product must be stored away from direct heat sources and sunlight, in a dry and clean place, under the conditions indicated on the label. **DISPOSAL:** the product does not contain hazardous substances and it may be disassembled and disposed of as regular municipal waste. **PRODUCT ORIGIN:** the last letter on the label of the product code REF indicates the source of the material (S = porcine; E = equine). Symbols used on the label: ☺ single use ☒ year and month of expiry [O] lot number [L] carefully read the instructions [E] Sterile. Sterilization performed by irradiation with gamma rays. If the packaging is open or damaged the sterility of the product is no longer guaranteed [T] limitation of the temperature

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