# LeadFIX®

**BIOVISION**  
BIOMATERIAL

BIOVISION GmbH  
Am Vogelherd 52  
98693 Ilmenau  
Germany  
info@biovision.de  
+49 3677 64 07 0  
+49 3677 64 07 13

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## EN INSTRUCTION MANUAL

### 1. Intended use

Fixing all kinds of membranes in dental surgery for a controlled regeneration of bones and tissues. It is recommended to check the suitability of the products for the planned intervention before use.

#### Indications

Fixation of membranes after e.g.:

- Fenestration defects at implant insertion
- Sinus floor elevation and augmentation
- Jaw ridge grafting
- Augmentation after block removal in the retromolar and chin region
- Augmentation of periodontal pockets

#### Contraindications

- Patients with pronounced osteoporosis
- Medications and therapies influencing bone wound healing
- Infections at the implantation site

### 2. Factors that can impair the success of surgery

- Insufficient bone quality (e.g. osteoporosis)
- Local bone tumors
- Systemic diseases and metabolic disorders
- Anamnestic indications of infections
- Drug addiction, alcohol and drug abuse
- Mental incapacity of the patient to understand and adhere to the instructions of the physician

#### Possible negative impacts

In rare cases, the following product- or procedure-related negative impacts can occur with the use of LeadFIX®:

- Foreign body reactions like local edema and osteolysis
- Infections
- Bleedings, Hematoma
- Wound dehiscence, fenestrations
- Swelling of the gums
- Heterotopic ossifications
- Dislocation and migration
- Allergic reaction
- Paresthesia
- retarded or partial healing of the wound

### 3. Patient information briefing

The physician is to inform the patient about possible impairments of the operation success and possible negative effects according to the sections "contraindications" and "possible negative effects", insofar as they are applicable for the patient. Furthermore, the patient is to be informed, which measures he or she can undertake to reduce possible impairments and effects.

### 4. Operation technique

A surgeon who is not familiar with the specific product is expected to study the surgery technique in detail before the operation. The product specific documentation is available from the manufacturer.

### 5. Handling

**NOTE:** This product is sterile packed in a protective packing. The LeadFIX is located inside the bottle. If you are unable to locate sufficient pins, please check the inside of the stopper as well.

**WARNING:** Danger of injury and risk of damage of products: The non-observance of these instructions for use can lead to injuries of patients, users and third parties, as well as to damage of the product. Carefully read through the relevant instructions for use and always adhere to the described instructions. Verify the function of the products that are used in combination.

**WARNING:** Danger of injury: Through the incorrect use of medical instruments there is the hazard of injury for patients. The users of medical instruments must have a corresponding medical qualification and be familiar with their application.

**WARNING:** Before using check the expiry date, as well as the packing for any damage. Products with elapsed expiry date, respectively damaged packing must not be used.

**WARNING:** Risk of infection: Through the use of non-sterile implants there is risk of infection for the patient. Only remove the implants from the sterile packing immediately before the application.

**WARNING:** Risk of infection: Reconditioning of the disposable articles can endanger the functionality and safety of the product and lead to infections. You must not recondition disposable articles. Dispose of disposable articles appropriately.

### 6. Handling of the accessories

**WARNING:** Risk of infection: These instruments are not supplied in sterile condition. Through the use of non-sterile instruments there is danger of infection for patients, users and third parties. Inspect instruments for visible contaminations. Visible contaminations indicate that the reconditioning has not been performed, or has not been performed correctly. Recondition the instruments using validated processes before the first use, as well as before and after every further use.

**WARNING:** Risk of infection: Not correctly reconditioned instruments cause patients, users and third parties to be at risk for infections, as well as the danger of functional impairments of the instrument. Adhere to the instructions on "Cleaning, disinfection, care and sterilisation of instruments" and the documents accompanying the product.

**NOTE:** The instructions on "Cleaning, disinfection, care and sterilisation of instruments" can be requested from the manufacturer.

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## EN INSTRUCTION MANUAL

**WARNING:** Danger of injury: An overload through too strong exertion of force can lead to fractures, deformation and functional impairments of the medical device and to injuries of the patient or user. Do not overload instruments. Do not bend deformed instruments back to the original position.

### 7. Package sizes and instruments

LeadFIX® is available in packages with 2 pins (PL 0121). The following instruments are available for LeadFIX®:

Article number	Instrument
PL 0114	Drill 0.8 mm
PL 0115	Drill 0.9 mm
PL 0116	Drilling template
PL 0119	Dispenser
PL 0120	Hammer
PL 0126	Applicator head straight
PL 0127	Applicator head angled
PL 0128	Applicator handle

### 8. Application

**ATTENTION:** The product must be exclusively implanted with the instruments intended for this purpose. During the implantation the patient must be treated and observed with the usual medical diligence. This includes the control of the progression of the procedure, and the monitoring of the vital parameters and the anaesthesia.

### 9. Traceability

The self-adhesive identification labels that are enclosed with the original packing must be stored for traceability and be adhered to the operation report, as well as the patient file. The individual at the hospital responsible for the stock keeping must maintain records of the available products, including article number (REF) and batch number (LOT), so that in the event of changes or recalls action can be quickly taken.

### 10. Packaging and sterility

The LeadFIX® are double bagged and sterilized by gamma rays with a dose of at least 25 kGy. The sterility is guaranteed until the date of expiry, as far as the packaging has neither been opened nor damaged.

### 11. Product material

The bioresorbable LeadFIX® is made of poly(L-co-D,L-lactate 70/30). This bioresorbable material is biocompatible and degradable.

### 12. Reusability

**WARNING:** Risk of infection: Reconditioning of disposable articles can endanger the functionality and safety of the product and lead to infections. You must not recondition disposable articles. Dispose of disposable articles appropriately.

### 13. Storage conditions

The product must be stored dry at max. 25 °C in the unopened original packing. Products with elapsed expiry date, respectively damaged packing must not be used.

### 14. Responsibility

Complications or other negative effects due to incorrect indication or surgical technique, unsuitable selection or handling of products or instruments, asepsis etc. fall under the responsibility of the surgeon.

### 15. Conformity with directives

This medical device is marked with CE labelling according to the Medical Device Directive (MDD) 93/42/EEC. If there is an identification number after the CE labelling then this designates the responsible specified authority.

### 16. Disposal

No special measures are necessary for the disposal of this medical device. The national laws and directives must be observed.

### 17. Packing symbols

	Manufacturer		The product is not suitable for reuse
	Observe the instructions for use		Do not resterilise
	Maximum shelf life of the product		Temperature limit
	Do not use the product with damaged packing		Store dry
	Sterilised through irradiation		

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