- (EN) see www.fkg.ch for additional information
 (FR) informations supplémentaires sur www.fkg.ch
 (DE) weitere Informationen auf www.fkg.ch


R-Motion®

Endodontic instruments
 Instruments endodontiques
 Endodontische Instrumente

FKG REF. 99.7AA.10.21A.XX - n°112 - 2021/10

R-Motion®

Description - Beschreibung

ISO diameter (1) identification

Identification diamètre ISO (1)

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(EN) General information

Parameters for use:

R-Motion® instruments must only be used in **reciprocating** (counterclockwise) mode.

Indications for use:

R-Motion® endodontic instruments are intended for use in medical or hospital facilities, by qualified health professionals.

R-Motion® Glider and R-Motion® Glider C instruments are designed for catheterisation of root canals in reciprocating mode.

R-Motion® 20/25/30/40/50 instruments are designed for shaping and cleaning root canals in reciprocating mode.

Contraindications:

Nickel-Titanium instruments must not be used on individuals with a known allergic sensitivity to these metals.

Adverse reactions: None known.

Notice and precautions:

- Take several radiographs from different angles to determine the anatomy of the root canals (length, width and curvature).
- Irrigate the canal thoroughly and frequently after the R-Motion® instrument is removed.
- Regularly inspect the R-Motion® instrument during use, and discard if it shows any signs of wear (for example : straightening).
- If the R-Motion® instrument is not progressing easily, withdraw it from the canal, clean it, check its cutting edges, then irrigate the canal and repeat the process with an ISO 010 hand file.
- Single use instruments must not be retreated for subsequent use. This could compromise their integrity.
- Always inspect the instrument(s) before use, and discard (it) them if there are any visible defect(s).
- When an instrument reaches the end of its life, please dispose of it in accordance with the applicable laws and regulations.

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1. Create straight-line coronal and radicular access.
2. Use an ISO 010 hand file to explore the canal.
3. Determine the working length using a radiograph and/or an apex locator.
4. Use the mechanised **R-Motion® Glider** instrument, applying gentle strokes of 2-3 mm to establish a glide path until the working length is reached. After 3 strokes, remove and clean the instrument, then irrigate; recapitulate with an ISO 010 hand file. *Use of the **R-Motion® Glider C** is recommended for canals which require an instrument with more rigidity, and/or if the **R-Motion® Glider** is no longer progressing and can not reach the working length.*
5. Select the R-Motion® single-file shaping instrument which is best suited to the anatomy of the canal being treated. In most cases, and depending on your approach to apical preparation, the **R-Motion® 25** or **R-Motion® 30** will be appropriate. *If the radiographs do not allow you to clearly identify the anatomy of the canal, it can be gauged using hand files.*
6. Start the shaping procedure, applying your usual irrigation/ disinfection protocol.
7. Use gentle inward-outward movements of 2-3 mm, applying very light apical pressure and allowing the instrument to progress passively along the canal. After 3 strokes, remove and clean the instrument, then irrigate; recapitulate with an ISO 010 hand file.
8. Repeat the previous step until the instrument has reached its working length. *If necessary, a brushing motion may be applied to remove coronal interferences and/or ensure uniform shaping of irregular canals. If the instrument is no longer progressing, make one or more passes with the **R-Motion® 20** instrument, then resume using the **R-Motion® 25** or the **R-Motion® 30**. Once the working length has been reached, remove the instrument to prevent over-enlargement of the foramen.*

Once all the canals have been shaped, proceed with the cleaning of the canals. *For an optimal result, the use of **XP-endo® Finisher** is recommended.*

Once the cleaning of the canals is completed, proceed to the next phase of treatment:

9. Obturate with gutta-percha and sealer. **TotalFill® BC Points™** and **TotalFill® BC Sealer™** are recommended.

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This is to certify that:

Emergo Global Consulting, LLC

2500 Bee Cave Road
Building 1, Suite 300
Austin
Texas
78746
USA

Holds Certificate Number:

FS 594789

and operates a Quality Management System which complies with the requirements of ISO 13485:2016 for the following scope:

Consulting services for regulatory and clinical affairs, quality system implementation and maintenance, training and related services, as well as in country representation services, such as European Authorized Representation and Australian Sponsor, delivered to the medical device industry.

For and on behalf of BSI:

Graeme Tunbridge, Senior Vice President Medical Devices

Original Registration Date: 2013-04-12

Latest Revision Date: 2023-03-15

Effective Date: 2023-04-02

Expiry Date: 2026-04-01

Page: 1 of 2

...making excellence a habit.™

Certificate No: FS 594789

Location	Registered Activities
Emergo Global Consulting, LLC 2500 Bee Cave Road Building 1, Suite 300 Austin Texas 78746 USA	Consulting services for regulatory and clinical affairs, quality system implementation and maintenance, training and related services, as well as in country representation services, such as European Authorized Representation and Australian Sponsor, delivered to the medical device industry. AMAR Registration and License Holder services delivered through Austin, Texas, USA and Matam, North Haifa, Israel locations.
I.L. Emergo Israel Ltd. Andrei Sakharov 9 Building 25 - 8th floor Matam North Haifa 3190501 Israel	AMAR Registration and License Holder services delivered through Austin, Texas, USA and Matam, North Haifa, Israel locations.
Emergo Europe Westervoortsedijk 60 6827 AT, Arnhem The Netherlands	Consulting services for regulatory and clinical affairs, quality system implementation and maintenance, training and related services, as well as in country representation services, such as European Authorized Representation delivered to the medical device industry.

Original Registration Date: 2013-04-12

Latest Revision Date: 2023-03-15

Effective Date: 2023-04-02

Expiry Date: 2026-04-01

Page: 2 of 2

This certificate was issued electronically and remains the property of BSI and is bound by the conditions of contract.

An electronic certificate can be authenticated [online](#).

Printed copies can be validated at www.bsigroup.com/ClientDirectory

Issuing Body: BSI Group The Netherlands B.V., John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands

BSI Group The Netherlands B.V. is registered in The Netherlands under number 33264284 | A Member of the BSI Group Holdings B.V.

Contact Office: 12950 Worldgate Drive, Suite 800, Herndon, VA 20170-6007 USA.