INTENDED USE


The device is meant to be used to maintain the access channel of the screw in a screw-retained prosthesis impervious so that it prevents the penetration of cement or filling material (no matter if provisional or final), thus avoiding bacterial colonization of free space and subsequent migration towards the tissues.

HOW TO USE

1. Wash your hands carefully before opening the blister.
2. Take the device out of the blister packaging using a dental tweezer.
3. Place the device in the hollow seat; then put SilverPlug® into the screw tunnel (Fig. 1).
4. Cut SilverPlug® in order to make it fit into the screw tunnel (Fig. 2).
5. Place SilverPlug® into the cavity and close with composite (Fig. 3).

WARNINGS

- Don't use SilverPlug® if the packaging is open or damaged.
- For single use only, please do not re-use the device in order to avoid possible contamination or lack of effectiveness.
- After opening the blister packaging, SilverPlug® must be used immediately.
- Keep SilverPlug® away from children.
- SilverPlug® is intended to be used only by trained medical personnel in suitable locations, namely dental laboratories and clinics.
- Dispose of used devices as special hazardous waste with infectious risk.

PRODUCT DESCRIPTION

SilverPlug® a preformed solid body. In the shape of an extended truncated cone easy to fit into the dental prosthetic housing seat of an Implant.

Once plugged in the seat, SilverPlug® reduces the free spaces that could be colonized by bacteria, thus reducing the eventual bacterial proliferation simply through the preservative action of the zeolite silver present in the device. SilverPlug® is obtained by using Thermo Plastic Elastomer (TPE) with silver in zeolite to the extent of 0.7% of the total weight. The silver contained in zeolite is between 1% and 3.0%.

The silver zeolite has the function to preserve the product given that it manages to maintain the microenvironment of the device, thus reducing the proliferation of bacteria. There is no release of silver zeolite into the body. The TPE resists the corrosive attack of both acidic and basic environment of the oral cavity and degrades the action of those pathogenic saprophytes and bacteria that might be found in the oral cavity. SilverPlug® replaces resins and other materials (such as gutta-percha, cotton, Teflon, wax, etc.) usually used to fill the cavity, but that generally don't guarantee the non-proliferation of bacteria and are even difficult to remove. SilverPlug® can be easily removed, in case it is necessary to act on the implant, with no risk of damaging the screw.

Medical device for occlusion in implant-prosthetic dentistry.

• Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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