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**en INSTRUCTIONS FOR USE**

**POLYAMIDE SUTURE MATERIALS**

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**MATERIALI DA SUTURA IN POLIAMMIDE**

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**MATERIALES DE SUTURA DE POLIAMIDA**

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- **NYLON**
- **SERALON®**
- **SUPRAMID**

**WARNINGS AND PRECAUTIONS**

- Read the instructions for use in their entirety before using the product.
- The product should only be used by qualified personnel, that is responsible to be familiar with both appropriate surgical techniques and anatomy prior to use.
- The product is delivered sterile (sterilization method: ethylene oxide; for SUPRAMID in cassette packs: gamma irradiation).

STERILE EO

STERILE R

- The product is intended for single use only.

This does not apply to cassette packs.

- For cassette packs, the shelf life after opening should be considered.
- Do not re-sterilize or re-use the product as critical changes in performance and safety are possible.

- Do not use after expiration date.

- Do not use if the product shows signs of damage (e.g. cuts or kinks) or if the outer packaging or sterile packaging has been damaged or unintentionally opened prior to use.

- Use the product only for the intended indication.

**DESCRIPTION AND INTENDED PURPOSE**

NYLON, SERALON®, and SUPRAMID are synthetic, non-absorbable, surgical suture materials made of polyamide. SUPRAMID in the cassette pack is supplied in an isopropanol-containing preserving solution.

Polyamide suture materials are available in various thread structures, gauge sizes and lengths. They are supplied in combinations with various types of needles made of surgical steel or without needles. The suture material complies with the current edition of the European Pharmacopoeia monograph (EP) "Sutures, sterile non-absorbable

(fila non resorbilia sterilia)". In deviation from this standard, the maximum thread length in the case of the cassette packs is not 4, but 100 m.

- **NYLON** black

NYLON is a monofilament thread and is dyed with Haematein (Logwood) black, Colour Index No. 75290. NYLON is available in very fine gauge sizes.

- **SERALON**<sup>®</sup> blue and undyed

SERALON<sup>®</sup> is a monofilament thread. SERALON<sup>®</sup> blue is dyed with Phthalocyanine blue, Colour Index No. 74160.

- **SUPRAMID** black and undyed

SUPRAMID is made of a multifilament twisted filament-bundle and a sheath. Therefore, it appears as a monofilament. SUPRAMID in gauge size EP1 and thinner is a monofilament. SUPRAMID black is dyed with pigment black, Colour Index No. 77266.

Polyamide threads are suitable in principle for all wound sutures and ligatures in which non-absorbable thread material is indicated and long-term stability of the suture material is not required. Threads left in tissues for a long time can lose tensile strength and eventually can fragment. NYLON can be manufactured to a minimum gauge size of EP 0.1 and is therefore suitable for sutures in microsurgery and ophthalmic surgery. Polyamide threads can also be used as holding sutures and for marking purposes.

## **INSTRUCTIONS FOR USE**

In the choice and use of suture material, account should be taken of the condition of the patient, the experience of the qualified personnel, the surgical technique, and the type and size of the wound. Standard knot-tying techniques should be used to ensure knot-holding security. In order to prevent needle damage during suturing it is recommended that the needle be grasped at an area between the end of its first third (as seen from its end) and its midpoint.

Unneeded suture material may be used for ligature or may be connected with surgical eye needles. For a list of suitable needles and application instructions, please refer to the Instructions for Use of Serag-Wiessner sterile eye needles. The suitability of other manufacturers' products has not been tested. When using other brands, their specific instructions must be observed.

## **ACTION**

The function of a medical suture in wound care is to attach tissue to tissue or tissue to foreign material. Ligatures are a special type of suture intended to close hollow organs.

## **CONTRAINDICATIONS**

The device is contraindicated for use in direct contact with the heart, the central circulatory system and the central nervous system. Because of the loss of tensile strength that occurs in the body, polyamide threads must not be used when prolonged stability of the suture material is required.

The device must not be used in patients with allergies to any of the product components.

## **RESTRICTIONS AND LIMITATIONS OF USE / NOTES / INTERACTIONS**

The requirements of the suture material and the techniques required vary with the indication. In the choice of suture material, account should be taken of the characteristics of the suture material concerned. In view of the high elasticity of thicker monofilament polyamide threads, it is recommended that one to two additional knots be applied in order to achieve knot-holding security. Like all foreign bodies, the device

can give rise to calculus formation when in pro-longed contact with salt-containing solutions (e.g. in the biliary or urinary tract).

When polyamide threads are used in ophthalmic surgery, the following warnings should be borne in mind: patients undergoing eye operations require prolonged observation; corneal polyamide threads should be removed after 3 to 12 months; laser irradiation accelerates the tissue breakdown of polyamide threads. Monofilament threads (SERALON<sup>®</sup>, NYLON) are preferable for use as sutures in infected tissue or tissue predestined for infection.

Surgical instruments should be handled with care in order to avoid damage to the threads. Bending of needles or grasping them outside of the recommended area can result in needle breakage. Because of the risk of infection from needle pricks, needles should be handled with particular care.

### **ADVERSE REACTIONS**

The following adverse reactions can occur in association with correct use of the device: local irritation, inflammatory reactions to the foreign body; formation of foreign-body giant cells, suture fistulas or granulomas.

Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

### **STORAGE CONDITIONS**

The storage conditions apply as indicated on the label.

### **INSTRUCTIONS FOR REMOVAL FROM THE PACK**

**Minibox:** For easy removal, the thread should be withdrawn from the Minibox not obliquely, but vertically upwards.

**Minibox - microsutures in protective sleeve:** Remove the sleeve from the minibox attachment; grasp the needle with the needle holder and withdraw the thread from the protective sleeve; if a **short thread** is required, it should be cut through together with the protective sleeve; in the case of **double-armed threads**, two individual threads can, if required, be withdrawn after the loop of thread at the end of the protective sleeve has been cut through.

**Cassette pack (also called suture dispenser or outer pack):** Remove the tamperproof closure (dispenser cap closure) and use the enclosed label to note the expiry date on the outer pack; observe the following instructions to guarantee aseptic removal for a period of 6 months after opening. Flip up the dispenser cap at the clip and fix it to the back part of the cap. If the outer pack has already been opened, i.e. the dispenser cap closure has been removed or opened, **disinfect** the inner surface of the dispenser **before use** with a commercially available alcohol-based product. Pull out the end of the suture with a sterile instrument. Discard the first 10 cm of the suture. Pull out the desired length of suture, cut it off flush, and close the cap. **After finishing the procedure, but no more than 4 hours later, or if required again disinfect** the inner surface of the dispenser and the cap with a commercially available alcohol-based product and close the cap.

### **DISPOSAL INFORMATION**

Disposal of unused products or waste material with particular care due to the potential risk of contamination and/or injury. In all other cases observe the specific national regulations.

## **MEANING OF THE SYMBOLS**

Symbols are explained at the end of the instructions for use.

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SERAG-WIESSNER GmbH & Co. KG