XP-endo® Rise Sequence

1 blister* includes / contient / enthält : 3 × (XP-endo® Rise Glider 15/04, XP-endo® Rise Shaper 30/04)

* Except for PromoSet XP-endo® Rise Sequence containing one sequence of two instruments.
 A l'exception des PromoSet XP-endo® Rise Sequence contenant qu'une séquence de deux instruments.
 Außer PromoSet XP-endo® Rise Sequence; dieses enthält eine Sequenz mit zwei Instrumenten.

- (EN)** see www.fkg.ch for additional information
- (FR)** informations supplémentaires sur www.fkg.ch
- (DE)** weitere Informationen auf www.fkg.ch

XP-endo® Rise Sequence

Endodontic instruments
 Instruments endodontiques
 Endodontische Instrumente

FKG REF. 99.X00.10.12A.XX - n°116 - 2021/12

XP-endo® Rise Sequence


Description - Beschreibung

XP-endo®
 Rise Glider
 15/04

Speed / vitesse / Drehzahl: 1000 rpm

Torque / couple / Drehmoment: 1 Ncm

XP-endo®
 Rise Shaper
 30/04

Speed / vitesse / Drehzahl: 1000 rpm

Torque / couple / Drehmoment: 1 Ncm

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(EN) General information

Indications for use:

These endodontic instruments are intended for use in medical or hospital facilities, by qualified health professionals. XP-endo® Rise instruments are designed for glide path and shaping of root canals.

Contraindications:

This product contains nickel and must not be used on individuals with a known allergic sensitivity to this metal.

Adverse reactions: None known.

Warnings and precautions:

- Take several radiographs from different angles to determine the anatomy of the root canals (length, width and curvature).
- Always inspect the instrument(s) before use, and discard (it) them if there are any visible defect(s).
- The number of canals to be treated with a XP-endo® Rise instrument is limited to a maximum of 8 root canals.
- Use the instruments XP-endo® Rise in the order of the sequence. Start with the glide path and then the shaping of the canal.
- Irrigate the canal thoroughly and frequently after the XP-endo® Rise instrument is extracted.
- Regularly inspect the XP-endo® Rise instrument during use, and discard if it shows any signs of wear (for example : straightening).
- If the XP-endo® Rise instrument is not progressing easily, remove it from the canal, clean it, inspect its cutting edges, then irrigate the canal and recapitulate with an ISO 010 hand file.
- When an instrument reaches the end of its life, please dispose of it in accordance with the applicable laws and regulations.
- These instruments have not been tested on children, pregnant nor breastfeeding women.
- Visit the FKG Dentaire SA website (www.fkg.ch) for all information relating to the applicable cleaning and sterilization methods.

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
Figure / Abbildung 2

Speed: 1000 rpm
Torque: 1 Ncm

1. Create straight-line coronal and radicular access.
2. Use an ISO 010 hand file to explore the canal.
3. Determine the working length using a radiograph. In addition, electronic working length determination can be achieved.
4. Use the mechanized XP-endo® Rise Glider instrument to perform glide path until the working length is reached.
The instrument XP-endo® Rise Glider are to be used with gentle strokes of 2-3 mm, applying very light apical pressure and allowing the instrument to progress passively along the canal. After 3 strokes, remove and clean the instrument, then irrigate; recapitulate with an ISO 010 hand file.
5. Start the shaping procedure, applying your usual irrigation protocol.
6. Insert the tip of the XP-endo® Rise Shaper into the canal until resistance (Fig. 1), retract (tip loose) and start the motor.
7. Use long gentle strokes to progress down to working length (WL) (Fig. 2). If WL is not reached in 8 strokes, stop, irrigate, and proceed again. Do not use pecking motion! Never force the instrument and always keep it spinning and moving while in the canal.
8. Once working length is reached, remove the instrument to avoid supplementary enlargement of the canal.
9. Irrigate the canal in order to eliminate suspended debris.
10. Choose a gutta percha point that is tight and seats at the correct WL. If 30/.04 is desired and does not fit, work the XP-endo® Rise Shaper for an additional 1-3 strokes.
11. Apply your usual disinfection protocol.

Once all the canals have been shaped, proceed with the deep cleaning (or disinfection) of the canals. *For an optimal result, the use of the Finisher is recommended.*

Once the cleaning of the canals is completed, proceed to the next phase of treatment :

12. Obturate with gutta percha and sealer. *TotalFill® BC Points™ and TotalFill® BC Sealer™ are recommended.*

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This is to certify that:

Emergo Global Consulting, LLC

2500 Bee Cave Road
Building 1, Suite 300
Austin
Texas
78746
USA

Holds Certificate Number:

FS 594789

and operates a Quality Management System which complies with the requirements of ISO 13485:2016 for the following scope:

Consulting services for regulatory and clinical affairs, quality system implementation and maintenance, training and related services, as well as in country representation services, such as European Authorized Representation and Australian Sponsor, delivered to the medical device industry.

For and on behalf of BSI:

Graeme Tunbridge, Senior Vice President Medical Devices

Original Registration Date: 2013-04-12

Latest Revision Date: 2023-03-15

Effective Date: 2023-04-02

Expiry Date: 2026-04-01

Page: 1 of 2

...making excellence a habit.™

Certificate No: FS 594789

Location	Registered Activities
Emergo Global Consulting, LLC 2500 Bee Cave Road Building 1, Suite 300 Austin Texas 78746 USA	Consulting services for regulatory and clinical affairs, quality system implementation and maintenance, training and related services, as well as in country representation services, such as European Authorized Representation and Australian Sponsor, delivered to the medical device industry. AMAR Registration and License Holder services delivered through Austin, Texas, USA and Matam, North Haifa, Israel locations.
I.L. Emergo Israel Ltd. Andrei Sakharov 9 Building 25 - 8th floor Matam North Haifa 3190501 Israel	AMAR Registration and License Holder services delivered through Austin, Texas, USA and Matam, North Haifa, Israel locations.
Emergo Europe Westervoortsedijk 60 6827 AT, Arnhem The Netherlands	Consulting services for regulatory and clinical affairs, quality system implementation and maintenance, training and related services, as well as in country representation services, such as European Authorized Representation delivered to the medical device industry.

Original Registration Date: 2013-04-12

Latest Revision Date: 2023-03-15

Effective Date: 2023-04-02

Expiry Date: 2026-04-01

Page: 2 of 2

This certificate was issued electronically and remains the property of BSI and is bound by the conditions of contract.

An electronic certificate can be authenticated [online](#).

Printed copies can be validated at www.bsigroup.com/ClientDirectory

Issuing Body: BSI Group The Netherlands B.V., John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands

BSI Group The Netherlands B.V. is registered in The Netherlands under number 33264284 | A Member of the BSI Group Holdings B.V.

Contact Office: 12950 Worldgate Drive, Suite 800, Herndon, VA 20170-6007 USA.