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**Bone and collagen matrix OsteoBiol® - English version -**

Class III medical device, bone and collagen matrix of heterologous origin, for the replacement of the autologous bone represented by the matrix and the porous form.

**Gen-Os®:** heterologous bone matrix and collagen, osteoconductive, resorbable, biocompatible, higroscopic. The product is moderately radiopaque. GRANULOMETRIA: 0,25-1 mm / 1-2 mm.

**CLINICAL INDICATIONS:** collagenated filler for non-infected, non-sclerotic, bone defects on non-load bearing and well-vascularized sites.

**PACKAGING:** sterile disposable glass vial or blister.

**DIRECTIONS FOR USE:** OsteoBiol® Gen-Os® must be hydrated prior to use with sterile solutions. For improved handling, OsteoBiol® Gen-Os® can be mixed with other devices such as OsteoBiol® TSV Gel. The mixture must be placed in the receiving site with sterile instruments.

**SIDE EFFECTS:** none known. Although the grafting procedure adopted has demonstrated the inhibition of adverse antibody reactions, the remote possibility of graft rejection cannot be excluded.

**CONTRADICTIONS:** do not use in patients who have pathologies or conditions that contraindicate surgery or who have an acute or chronic infection (osteomyelitis). Not indicated for patients with immunodeficiency or suffering from autoimmune diseases, particularly children and the elderly. Do not use directly or near large vessels or nerves. Do not use on poorly vascularized, necrotic or infected receiving sites. Do not use on sites that do not guarantee the stability and protection of the product after grafting. The product must not be used if the packaging is damaged or opened before use.

**WARNING:** The product has not been tested on pregnant women. The product is for single use only and, as such, should not be sterilized. The graft site must be properly washed and cleaned, and the product must be well stabilized and protected.

Product reserved for surgeons with experience in biomaterial implants: strictly adhere to the instructions given in the scientific literature on heterologous grafts, especially the scientific literature on OsteoBiol® devices available on the website [www.osteobiol.com](http://www.osteobiol.com), as well as the instructions herein. For every single application, the specialist surgeon has the responsibility to ensure proper clinical management of the biomaterials and the patient.

**EXPIRY DATE:** the product is valid for five years from the date of sterilization and the expiry date is indicated on the label.

**DISPOSAL:** the product does not contain hazardous substances and it may be disassembled and disposed of as regular municipal waste.

**PRODUCT ORIGIN:** the last letter on the label of the product code REF indicates the source of the material (S=porcino; E=equino).

**Symbols used on the label:** ☐ single use ☐ year and month of expiry ☐ carefully read the instructions ☐ Sterile. Sterilization performed by irradiation with gamma rays. If the packaging is open or damaged the sterility of the product is no longer guaranteed ☐ limitation of the temperature

**STOP/AGE:** The product must be stored away from direct heat sources and sunlight, in a dry and clean place, under the conditions indicated on the label.

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